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Application of capillary isoelectric focusing with whole column imaging detection in protein analysis

Li Bin, Ouyang weimin, Yu Shijian

(Beijing Lumiere Tech Limited, Beijing, 100080, China)

Abstract: Capillary isoelectric focusing - Whole column imaging detection is a protein detection technology developed in 1992, It is a strong powerful tool for protein drug analysis and food protein research based on its fast analysis speed, good reproducibility and accurate determination.

Key words: capillary isoelectric focusing; whole column imaging detection; Protein drugs; Food protein

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- [1] HJERTENS, ZHUMD. Adaptation of the Equipment for High Performance Electrophoresis to Isoelectric Focusing [J]. *Journal of Chromatography*, 1985, 346: 265—270.
- [2] 吴家齐. 毛细管等电聚焦电泳(cIEF)及其理论模拟计算新进展 [J]. *计算机应用于化学*, 2007, 24(1): 11-15.
- [3] WUJQ, PAWLISZYNJ. High Performance Capillary Isoelectric Focusing with a Concentration Gradient Detector [J]. *Analytical Chemistry*, 1992, 64: 219—224.
- [4] WUJQ, PAWLISZYNJ. Capillary Isoelectric Focusing with a Universal Concentration Gradient Imaging System Using a Charge Coupled Photodiode Array [J]. *Analytical Chemistry*, 1992, 64: 2934—2941.
- [5] WUJQ, PAWLISZYNJ. Dual Detection for Capillary Isoelectric Focusing with Refractive Index Gradient and Absorption Imaging Detectors [J]. *Analytical Chemistry*, 1994, 66: 867—873.
- [6] WUJQ, PAWLISZYNJ. Diode Laser Based Concentration Gradient Imaging Detector for Capillary Isoelectric Focusing [J]. *Analytica Chimica Acta*, 1995, 299: 337—342.
- [7] WUJQ, PAWLISZYNJ. Absorption Spectra and Multi capillary Imaging Detection for Capillary Isoelectric Focusing Using a Charge Coupled Device Camera [J]. *Analyst*, 1995, 120: 1567—157.
- [8] WUXZ, WUJNQ, PAWLISZYNJ. Fluorescence Imaging Detection for Capillary Isoelectric Focusing [J]. *Electrophoresis*, 1995(16): 1474—1478.
- [9] MAOQL, PAWLISZYNJ. Capillary isoelectric focusing with whole column imaging detection for analysis of protein and peptides [J]. *Journal of Biochemical and Biophysical Methods*, 1999, 39: 93—110.