



Designation: C1001 – 19

Standard Test Method for Radiochemical Determination of Plutonium in Soil by Alpha Spectroscopy¹

This standard is issued under the fixed designation C1001; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon (ϵ) indicates an editorial change since the last revision or reapproval.

1. Scope

1.1 This test method covers the determination of plutonium in soils at levels of detection dependent on count time, sample size, detector, background, and tracer yield. This test method describes one acceptable approach to the determination of plutonium in soil.

1.2 This test method is designed for 10 g of soil, previously collected and treated as described in Practices C998 and C999, but sample sizes up to 50 g may be analyzed by this test method. This test method may not be able to completely dissolve all forms of plutonium in the soil matrix.

1.3 The values stated in SI units are to be regarded as standard. The values given in parentheses are for information only.

1.4 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety, health, and environmental practices and determine the applicability of regulatory limitations prior to use.* Specific precautionary statements are given in Section 10.

1.5 *This international standard was developed in accordance with internationally recognized principles on standardization established in the Decision on Principles for the Development of International Standards, Guides and Recommendations issued by the World Trade Organization Technical Barriers to Trade (TBT) Committee.*

2. Referenced Documents

2.1 *ASTM Standards:*²

C859 Terminology Relating to Nuclear Materials

C998 Practice for Sampling Surface Soil for Radionuclides

C999 Practice for Soil Sample Preparation for the Determination of Radionuclides

C1163 Practice for Mounting Actinides for Alpha Spectrometry Using Neodymium Fluoride

C1284 Practice for Electrodeposition of the Actinides for Alpha Spectrometry

D1193 Specification for Reagent Water

D3084 Practice for Alpha-Particle Spectrometry of Water

D3648 Practices for the Measurement of Radioactivity

D7282 Practice for Set-up, Calibration, and Quality Control of Instruments Used for Radioactivity Measurements

3. Terminology

3.1 Except as otherwise defined herein, definition of terms are as given in Terminology C859.

4. Summary of Test Method

4.1 Plutonium is extracted from the soil with a mixture of nitric, hydrofluoric, and hydrochloric acids in the presence of ²⁴²Pu or ²³⁶Pu isotopic tracer (See Appendix for purification and standardization of ²³⁶Pu tracer). Plutonium is isolated by anion exchange, then electrodeposited onto a polished metal disk for determination by alpha spectrometry. As an option, the plutonium may be prepared for alpha spectrometry measurement by using coprecipitation with neodymium fluoride. The range of chemical yield is between 40 and 90 %. The test method is based on a published procedure (1).³

5. Significance and Use

5.1 A soil sampling and analysis program provides a direct means of determining the concentration and distribution of radionuclides in soil. A soil analysis program has the most significance for the preoperational monitoring program to establish baseline concentrations prior to the operation of a nuclear facility. Soil analysis, although useful in special cases involving unexpected releases, is a poor technique for assessing small incremental releases and is therefore not recommended as a method for routine monitoring releases of radioactive material. Nevertheless, because soil is an integrator

¹ This test method is under the jurisdiction of ASTM Committee C26 on Nuclear Fuel Cycle and is the direct responsibility of Subcommittee C26.05 on Methods of Test.

Current edition approved Nov. 1, 2019. Published December 2019. Originally approved in 1983. Last previous edition approved in 2011 as C1001 – 11. DOI: 10.1520/C1001-19.

² For referenced ASTM standards, visit the ASTM website, www.astm.org, or contact ASTM Customer Service at service@astm.org. For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

³ The boldface numbers in parentheses refer to the list of references at the end of this standard.

and a reservoir of long-lived radionuclides, and serves as an intermediary in several of the plutonium pathways of potential importance to humans, knowledge of the concentration of plutonium in soil is essential.

6. Apparatus

6.1 *Electrodeposition Apparatus (2)*, see Practice [C1284](#).

6.2 *Alpha Spectrometer*, capable of 40 to 50 keV resolution on actual samples electrodeposited on a flat, mirror-finished metal planchet, and a fractional counting efficiency greater than 17 %, and a background less than 0.010 cpm over each designated energy region. Resolution is defined as the full width half maximum (FWHM) in keV, the distance between those points on either side of the alpha peak where the count is equal to one-half the maximum count. Also see Practices [D3084](#) and [D3648](#).

NOTE 1—A regular program of measurement control operations should be conducted for the alpha spectrometry system, such as regular background checks, source check to determine system stability, control charting, and careful handling of samples during changing. See Practice [D7282](#).

6.3 *PTFE-polytetrafluoroethylene beakers*, 250-mL.

6.4 *TFE-fluorocarbon stirring rod*.

6.5 *Centrifuge bottle*.

6.6 *Centrifuge*.

6.7 *Hot water bath*.

6.8 *Anion exchange columns*.

7. Reagents

7.1 *Purity of Reagents*—Reagent grade chemicals shall be used in all tests. Unless otherwise indicated, it is intended that all reagents shall conform to the specifications of the Committee on Analytical Reagents of the American Chemical Society, where such specifications are available ([3](#)). Other grades may be used, provided it is first ascertained that the reagent is of sufficiently high purity to permit its use without lessening the accuracy of the determination.

7.2 *Purity of Water*—Unless otherwise indicated, references to water shall be understood to mean reagent water as defined in Specification [D1193](#), Type III.

7.3 Reagent blanks should be run to determine their contribution to the sample result.

7.4 *Ammonium Hydroxide (sp gr 0.90)*—Concentrated ammonium hydroxide (NH_4OH).

7.5 *Ammonium Iodide (NH_4I) (1 M)*—Dissolve 14.5 g of NH_4I in 100 mL water.

7.6 *Ammonium Iodide, Hydrochloric Acid Solution ($\text{NH}_4\text{I-HCl}$)*—Add 25 mL 1 M ammonium iodide to 500 mL concentrated hydrochloric acid. Prepare fresh prior to use.

7.7 *Analytical Grade Anion Exchange Resin Type 1, 8 % cross-linked, 100 to 200 mesh, chloride form*—Store the resin in water in a wash bottle, transfer sufficient resin to a 1.3-cm ID ion exchange column to give a 10-cm bed of settled resin.

Convert the resin to the nitrate form by passing 100 mL of HNO_3 (8M) through the column at maximum flow rate.

7.8 *Boric Acid (H_3BO_3)*.

7.9 *Hydrochloric Acid (sp gr 1.19)*—Concentrated hydrochloric acid (HCl).

7.10 *Hydrochloric Acid (3M)*—Add 250 mL of concentrated HCl (sp gr 1.19) to 500 mL of water, mix, and dilute to 1.0 L with water.

7.11 *Hydrochloric Acid (1.7M)*—Add 142 mL of concentrated HCl (sp gr 1.19) to 500 mL of water, mix, and dilute to 1.0 L with water.

7.12 *Hydrofluoric Acid (48 to 51 %)*—Concentrated hydrofluoric acid (HF).

7.13 *Iron Carrier solution (10 g Iron (III)/L)*—Dissolve 10.0 g iron metal in HCl (1.7M) and dilute to 1 L with HCl (1.7M).

7.14 *Nitric Acid (sp gr 1.42)*—Concentrated nitric acid (HNO_3).

7.15 *Nitric Acid (8M)*—Add 500 mL of concentrated HNO_3 (sp gr 1.42) to 400 mL of water, mix, and dilute to 1.0 L with water.

7.16 *Nitric Acid (1.8M)*—Add 112 mL of concentrated HNO_3 (sp gr 1.42) to 500 mL of water, mix, and dilute to 1.0 L with water.

7.17 *Octyl alcohol*.

7.18 *Sodium Bisulfite (NaHSO_3)*.

7.19 *Sodium Hydroxide (50 %)*—Dissolve 500 g of NaOH in 500 mL water. A 50 % NaOH solution is available commercially.

7.20 *Sodium Nitrite (NaNO_2)*.

7.21 *Thymol Blue Indicator, Sodium Salt, 0.02 % Solution*.

7.22 *National Standard Traceable ^{236}Pu Reagent*.

7.23 *National Standard Traceable ^{242}Pu Reagent*.

8. Sampling

8.1 Collect the sample in accordance with Practice [C998](#).

8.2 Prepare the sample for analysis in accordance with Practice [C999](#).

8.3 Samples consisting of 10 to 50 g of soil can be readily analyzed by the procedure. In order to obtain more representative samples, as well as lowering the minimum detectable concentration, the analyses of large soil samples is desirable. In general, it is poor practice to use less than 10 g of sample, unless replicate analyses are performed, because of needed sensitivity to determine lower levels of activity.

9. Calibration and Standardization

9.1 The fractional counting efficiency of the alpha spectrometer is used to determine the chemical yield. The fractional efficiency of the alpha spectrometer is determined as the ratio of the observed count rate to the known disintegration rate.

9.2 The quantity of the tracer to be used should be in the expected range (but not less than 0.17 Bq) of the isotopic

activity being determined so that the statistical uncertainty in the chemical yield determination will not be larger than that of the nuclide being determined.

10. Precautions

10.1 Adequate laboratory facilities, such as fume hoods and controlled ventilation, along with safe techniques, must be used in this procedure. Extreme care should be exercised in using hydrofluoric and other hot, concentrated acids. Use of proper gloves is recommended. Refer to the laboratory's chemical hygiene plan and other applicable guidance for handling chemical and radioactive materials and for the management of radioactive, mixed, and hazardous waste.

10.2 Hydrofluoric acid is a highly corrosive acid that can severely burn skin, eyes, and mucous membranes. Hydrofluoric acid differs from other acids because the fluoride ion readily penetrates the skin, causing destruction of deep tissue layers. Unlike other acids that are rapidly neutralized, hydrofluoric acid reactions with tissue may continue for days if left untreated. Familiarization and compliance with the Safety Data Sheet is essential.

11. Procedure

11.1 Weigh a 10.0 ± 0.1 g soil aliquot into a 250 mL PTFE beaker.

11.2 Wet sample with distilled water and add a known quantity of ^{236}Pu or ^{242}Pu tracer.

11.3 Add concentrated HNO_3 (sp gr 1.42) a few drops at a time as fast as the frothing and vigor of the reaction will permit until the entire sample is covered.

11.4 Add 60 mL more of the concentrated HNO_3 (sp gr 1.42) and 30 mL of concentrated HF (48 to 51 %) and digest on a hotplate with frequent stirring (TFE-fluorocarbon stirring rod) for about 1 h (Note 2 and Note 3).

NOTE 2—For organic soils, first add the nitric acid only in small portions while stirring. If the solution threatens to overflow as a result of froth generation, add a few drops of octyl alcohol and stir. Digest on a hotplate until the evolution of reddish-brown fumes is reduced to a barely visible level. Cool to room temperature before carefully adding the concentrated hydrofluoric acid (48 to 51 %) and digesting for an hour.

NOTE 3—For larger soil aliquots, larger amounts of the acids (in the same proportions) should be used. For example, for a 50 g sample, use 200 mL concentrated HNO_3 and 100 mL HF, etc., with appropriately sized containers.

11.5 Remove from the hotplate and cool somewhat before adding 30 mL concentrated HNO_3 (sp gr 1.42) and 30 mL concentrated HF (48 to 51 %). Digest on the hotplate with intermittent stirring for an additional 1 h.

11.6 Remove from the hotplate and cool. Carefully add 20 mL concentrated HCl (sp gr 1.19) and stir. Heat on hotplate for 45 min with occasional stirring.

11.7 Add about 5 g of powdered boric acid and digest for an additional 15 min with occasional stirring.

11.8 Add approximately 200 mg of sodium bisulfite and continue heating until the solution has evaporated to a liquid volume of approximately 20 mL.

11.9 Add 50 mL of water and digest on a hotplate while stirring for 10 min to dissolve soluble salts.

11.10 Cool and transfer approximately equal parts of the total sample into centrifuge bottles with a minimum of water from a wash bottle. If equipment for large volume centrifugation is not available, the two precipitations in 11.11 – 11.21 may be performed in a beaker, allowing the precipitate to settle, decanting the supernate, and then completing the separation by centrifugation on a smaller scale.

11.11 Add 1.0 mL of iron carrier solution (10 mg Fe (III)/mL) to each centrifuge bottle and mix (Note 4).

NOTE 4—It may not be necessary to add the iron carrier if a sufficient amount of iron is present in the soil.

11.12 Add NaOH (50 %) with mixing to each bottle to a pH of about 9 (using pH paper) to precipitate ferric hydroxide. Add 5 to 10 mL excess NaOH and mix for 1 min.

11.13 Centrifuge for approximately 5 min, decant, and discard the supernate.

11.14 Dissolve each precipitate with about 30 mL HNO_3 (8M) (60 mL total) saturated with boric acid. (Approximately 7 g of boric acid/30 mL HNO_3 (8M).) Digest in a hot water bath for 10 min.

11.15 Cool and centrifuge for approximately 5 min. Decant the supernate into the original 250 mL PTFE beaker and save.

11.16 Wash each residue with approximately 10 to 20 mL (20 to 40 mL total) of HNO_3 (8M) saturated with boric acid. Centrifuge for 5 min and combine the supernates with that in 11.15. Discard the undigested soil.

11.17 Heat the supernate on a hotplate and evaporate to approximately 5 mL.

11.18 Add approximately 30 mL water and heat to dissolve the salts. Cool and transfer into a centrifuge tube.

11.19 Add concentrated ammonium hydroxide dropwise with mixing to a pH of approximately 9 (using pH paper) to precipitate ferric hydroxide.

11.20 Centrifuge and discard the supernate.

11.21 Dissolve the precipitate with a volume of concentrated nitric acid approximately equal to the volume of the precipitate and transfer using nitric acid (8M) into a 250-mL beaker. Add nitric acid (8M) to a total volume of approximately 75 mL. If the volume of the hydroxide precipitate is considerably greater than should be expected from the 10 mg of Fe(III) added, the final volume should be brought up to approximately 100 mL with nitric acid (8M) or, alternatively, the dissolved hydroxides should be evaporated to salts before the addition of the nitric acid (8M) solution. The final nitric acid molarity of the solution is not extremely critical, but should be in the range of 7 to 9.

11.22 Add approximately 200 mg of sodium nitrite crystals and stir with a stirring rod. Bring to a quick gentle boil on a hotplate and cool. Avoid prolonged heating.

11.23 Prepare an anion-exchange column as described in 7.7.