



## Standard Practice for Vacuum Distillation of Solvents From Solvent-Reducible Paints For Analysis<sup>1</sup>

This standard is issued under the fixed designation D 3272; the number immediately following the designation indicates the year of original adoption or, in the case of revision, the year of last revision. A number in parentheses indicates the year of last reapproval. A superscript epsilon ( $\epsilon$ ) indicates an editorial change since the last revision or reapproval.

*This standard has been approved for use by agencies of the Department of Defense.*

### 1. Scope

1.1 This practice covers a vacuum distillation procedure that separates the solvents from the nonvolatile portion of paints so they may be analyzed.

1.2 An alternative procedure for the analysis of solvents in paints is described in Practice D 3271.

1.3 The values stated in SI units are to be regarded as the standard. The values given in parentheses are for information only.

1.4 *This standard does not purport to address all of the safety concerns, if any, associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.*

### 2. Referenced Documents

#### 2.1 ASTM Standards:<sup>2</sup>

D 3271 Practice for Direct Injection of Solvent-Reducible Paints into a Gas Chromatograph for Solvent Analysis

### 3. Summary of Practice

3.1 The sample of whole paint with added tricresyl phosphate is heated by means of a silicone oil bath. The distillation is allowed to proceed for a period of time and then vacuum is applied. Finally dry air is passed through the system to facilitate the removal of any high-boiling solvents present.

### 4. Apparatus

4.1 Apparatus<sup>3</sup> is shown in Fig. 1. It is essential in assembling this apparatus, that the volume of the portion preceding the flask be kept to a minimum by using diameters of

connections and distances as small as possible. Include a gage in the assembly for measuring pressure.

4.2 *Vacuum Pump*, capable of pulling 2 mm Hg or less of vacuum pressure.

4.3 *Septum Sample Vial*, PTFE-fluorocarbon-faced.<sup>4</sup>

4.4 *Safety Shield*.

### 5. Reagents

5.1 *Methylene Chloride or Acetone*.

5.2 *Silicone Oil*.<sup>5</sup>

5.3 *Tricresyl Phosphate*.

### 6. Procedure

6.1 Immerse the receiver in a bath containing dry ice with methylene chloride or acetone of sufficient depth that the liquid reaches to the bottom level of the rubber stopper.

6.2 Pour approximately 15 g of the whole paint sample into the 50-mL distillation flask. Immediately follow with 10 mL of tricresyl phosphate and without delay mix the contents of the flask.

6.3 Assemble the apparatus and close the needle valve and stopcock attached to the capillary tubing. Immerse the distillation flask by raising a silicone oil bath preheated to 155 to 160°C (**Warning**—Care should be taken to prevent the temperature of the silicone oil exceeding 160°C because of possible formation of explosive mixtures of air and solvents from the sample being analyzed.) until it reaches the side arm (Note 1). Allow the distillation to proceed for 15 min without vacuum.

NOTE 1—The liquid level in the oil bath should be as high as is practical.

6.4 After 15 min apply vacuum slowly and carefully to avoid bumping or rapid distillation until 2 mm Hg is reached.

<sup>1</sup> This practice is under the jurisdiction of ASTM Committee D01 on Paint and Related Coatings, Materials, and Applications and is the direct responsibility of Subcommittee D01.21 on Chemical Analysis of Paints and Paint Material.

Current edition approved Dec. 1, 2003. Published December 2003. Originally approved in 1973. Last previous edition approved in 1976 as D 3272 – 76 (1998).

<sup>2</sup> For referenced ASTM standards, visit the ASTM website, [www.astm.org](http://www.astm.org), or contact ASTM Customer Service at [service@astm.org](mailto:service@astm.org). For *Annual Book of ASTM Standards* volume information, refer to the standard's Document Summary page on the ASTM website.

<sup>3</sup> Obtainable from SGA Scientific, Inc., 735 Broad St., Bloomfield, NJ 07005.

<sup>4</sup> Septum vials available from Precision Sampling Corp., Baton Rouge, LA 70315, have been found satisfactory for this purpose.

<sup>5</sup> The sole source of supply of silicone oil 710 known to the committee at this time is Dow Corning Corp., Midland, MI. If you are aware of alternative suppliers, please provide this information to ASTM International Headquarters. Your comments will receive careful consideration at a meeting of the responsible technical committee,<sup>1</sup> which you may attend.